

# BARTENDER™

BY: REVHECHT CONCEPTS LLC.

"Raising the Bar in Your Gym with The BarTender™"  
"The Original Bar Cane"

The BarTender™ is an industry specific tool designed to assist in the setting changes of the high rail on the women's gymnastics uneven bar apparatus. It is also used as a "ring still" tool for stilling the rings on the mens rings apparatus. It is designed to be used in both training and competitive settings.



## PROCEDURE:

1. Loosen all four cable tensioners of the apparatus.
2. Position the hook end of the BarTender™ onto the metal collar at the end of the high bar rail.  
**WARNING: To prevent marring and scratching of the high bar rail, only adjust the rail by the metal collar. RevHecht Concepts, LLC., will not be liable for injury or damage caused by failure to follow warnings and improper use of the BarTender™.**
3. Loosen the "T" handle lock.
4. Pull out the spring lock pin to release the high rail upright post. (Be sure you prepare to support the weight of the high rail assembly with the BarTender™).
5. Raise or lower the high rail to your desired setting.
6. Release the spring lock pin and assure that it is fully engaged into the appropriate notch on the high bar rail mounting upright post.
7. Tighten the "T" handle lock.
8. Tighten all four cable tensioners.
9. The bar cane can be stored on the apparatus by hanging the cane end from the "T" handle lock on the outside of the apparatus spreaders where it can be easily reached.

## SAFETY:

"Free Dropping" is the practice of quick, intermittent releasing of the spring lock pin to lower the high rail assembly. The weight of the high rail assembly falling abruptly onto the internal spring lock pin can create excessive wear and damage to the pin which can increase the risk of equipment failure. If used properly, the BarTender™ eliminates free dropping.

No need to grasp the notched inner upright when raising or lowering the bar which can pinch fingers and hands.

## CONVENIENCE:

No need to stand on mats, blocks, springboards, stools, ladders or chairs to adjust the high rail. (no need to climb up and stand on the spreader cross bar). Collapses to 26" and extends to 40". The BarTender™ can be conveniently hung on the apparatus for quick access.

## AVAILABILITY & PRICING:

\$136.00 + Tax or 10% off of two or more units.  
Shipping available (Prices subject to change).  
1 Year warranty on materials and workmanship.

## CONTACT INFO:

David Sedgwick - Butler Gymnastics Club, Butler, PA 16001  
Email: skillcoach1@yahoo.com • Phone: 412-860-9236  
Click the QR code or visit the website: [revhechtconcepts.com](http://revhechtconcepts.com)

The **BARTENDER™**  
Raising the bar in your gym

